Microderm Treatment

Weekly Skin Renewal

Device Instructions

**PLEASE READ THIS INSTRUCTION BOOKLET
CAREFULLY BEFORE USING THIS DEVICE.**

WARNINGS

- Do not use this device without first viewing the instructional video. Failure to use this product strictly in accordance with the instructions and the instructional video could result in injury to your skin.
- Do not operate this device without the plastic safety cap securely in place.
- Using the device in one place for too long may cause skin irritation and abrasions.
- Never use the device above the cheekbone or below the eyebrow.
- Avoid areas with broken capillaries, rosacea, scabbing or open sores, tender acne lesions, or sensitive areas of the skin.
- Not recommended for severe acne.
- Do not use any other type of mechanical exfoliation or chemical exfoliation (harsh acids, peels, etc.) for at least 4 days after using the Defy Microderm Treatment.
- Do not use products containing retinol or Retin-A for at least 4 days after using the device.
- Do not allow water to get inside the vacuum hole.
- It is best to consult a physician prior to use if you are unsure or have any questions about microderm use.

Defy Microderm Treatment Device Instructions

STEP 1: ASSEMBLY

There are four discs available. Use the white disc the first 2-3 times you use the device and move up as you become more comfortable using it. Those with sensitive skin should not use the black or red intensity levels.

Disc Color	Level of Intensity
White	Mild
Grey	Moderate
Black	Coarse
Red	Very Coarse

Important: There are two disc sizes—the larger for the body (i.e. elbows, hands, décolleté, etc.) and the smaller for the face. The large disc should not be used on the face.

1. Push disc tightly onto the metal tip at the end of the device. (To remove a disc, just pull off.)
2. Place the properly-sized cap over the disc. Screw until tight. The seal must be tight to ensure proper suction from the vacuum.
3. Plug cord into the device and twist $\frac{1}{4}$ turn to secure. Plug the cord into a power outlet.

Remember not to operate this device without the plastic safety cap securely in place. The top of the disc should be just barely below the top of the cap so that the vacuum can suck the skin to the disc.

STEP 2: PREPARATION

1. Cleanse and dry the desired treatment area, removing all dirt, makeup, and oil. Skin must be dry, not wet or damp, when using the device.
2. Make sure hair is pulled back and all facial accessories are removed.

STEP 3: TREATMENT

1. Using the fingers of your free hand, stretch the skin to create a smooth surface where the tip will make contact.
2. Turn the device on. Ensure that the cap is flush with the skin as you make contact. The motor sound will change when suction occurs.
3. Move the disc in upward and horizontal motions until you have treated all of the desired area.
4. Keep the tip moving constantly over skin surface—do not hover in one spot.
5. Use light pressure with early treatments, starting with one pass of

the device across the skin. As you become more familiar with the device, you can increase the pressure and/or number of passes according to your needs. Be careful not to apply too much pressure though. Flaky dead skin may appear on the surface of the skin you treat. Pinkness of the skin is normal because of the professional quality of this treatment.

STEP 4: DEVICE CLEAN UP

1. Rinse the filter (which keeps the dead skin cells from being sucked into the vacuum) thoroughly under warm running water.
2. Do not get water in the vacuum. It is important to clean this filter every time you use the device—if you don't, dead skin cells will clog the vacuum and the device will stop working properly.
3. Let the filter dry completely before putting the device away.
4. Use a soft cloth or tissue to wipe any residue out of the cap and motor shaft.

STEP 5: SKIN CLEANSE

Wipe or rinse any dead skin off your face. Apply Defy Triple Action Toner to better prepare skin for Defy treatment products.

STEP 6: FACIAL CARE PRODUCTS

Apply your other Defy facial care products as usual. Be sure to wear sunscreen during the day and avoid excessive sun exposure for 3-4 days following the treatment.

STEP 7: MAINTENANCE

The discs should be replaced every 3-4 treatments and are available for purchase at morinda.com/defy. For maximum results, grit on the discs must be fresh. The filter (available at morinda.com/defy) should also be replaced every 6 months.

Use the Defy Microderm Treatment device once a week for 10-12 weeks. After that, only use every other week for 2-3 months. The Defy Microderm Treatment can then be used once a month to maintain results.

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WARRANTY AND GUARANTEE:

The Defy Microderm Treatment comes with a 90-day warranty and 30-day money back guarantee. If you are in any way dissatisfied with your purchase, please send back the Defy Microderm Treatment and all of the contents, including device instructions and DVD, and we will refund you the total cost of the purchase (minus the cost of shipping).

The 90-day warranty period begins on the date the device was purchased. We cannot honor the warranty or guarantee if the device was damaged due to tampering, neglect, or misuse.

Questions? Call 1-800-445-2969.

Morinda, Inc.

333 West River Park Drive, Provo UT 84604

morinda.com/defy

This manual should not be copied, modified or distributed without written permission from Tahitian Noni International.

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